

# **NANOVEA** **OPTICAL** **PROFILERS**

---

*The New Standard  
of Profilometry*



# **NANOVEA<sup>®</sup>**

A Better Measure.

Offering more than **25 Years** of Material Science Experience



## **RESEARCH AND CONSULTATION**

Extensive range of research content such as brochures, application notes, publications, and videos.



## **EXPERT ASSISTANCE**

Dedicated Profilometry experts happy to guide you through any question or project request.



## **CUTTING EDGE INNOVATION**

At Nanovea we are always developing cutting edge technologies and standards. We innovate our instruments so that you can innovate your own products.



## **PRE AND POST INSTALLATION SUPPORT**

Full walk-through and guide to make sure the instrument is installed perfectly. Dedicated support team to help you after your instrument has been installed.

# ***INSTRUMENTS***



# ST400 OPTICAL PROFILER

- 200 x 150mm XY stages
- Video imaging integration
- Ideal for wide range of samples with varied geometries
- Chromatic confocal sensors w/ speed up to 200 times faster
- Rotational stage parallel or perpendicular to the testing plate
- Height sample clearance up to 200mm



## 3D CHROMATIC SENSORS



Standard Sensor



High Speed Sensor

## VIDEO OPTIONS




Atomic Force Microscope




Zoom Microscope

# THE STANDARD FOR PROFILOMETRY



**X - Y SCAN AREA**  
200 x 150mm Motorized



**HEIGHT RANGE**  
2.5nm to 25mm



**DESKTOP DIMENSIONS**  
62 x 62 x 82cm



**SCAN SPEED**  
40mm/s

# ST500 LARGE AREA OPTICAL PROFILER

- High speed large area measurement w/ high speed sensor
- 400 mm XY axis travel with a maximum speed up to 200 mm/s
- Video zoom camera to provide automated functions
- Measurements with a user friendly desktop platform



## 3D CHROMATIC SENSORS



Standard Sensor



High Speed Sensor

## VIDEO OPTIONS





Atomic Force Microscope




Zoom Microscope

# HIGH SPEED AND LARGE AREA MEASUREMENT



**X - Y SCAN AREA**  
400 x 400mm Motorized



**HEIGHT RANGE**  
2.5nm to 25mm



**DESKTOP DIMENSIONS**  
97 x 72 x 92cm



**SCAN SPEED**  
200mm/s

# JR25 PORTABLE OPTICAL PROFILER

- First truly portable non contact profilometer
- Weight less than 5.5 kg
- Lab quality results on the go
- Measurement capabilities up to 25mm x 25mm
- Able to measure samples at difficult angles
- Possible integration into automated robot arms and other equipment



## 3D CHROMATIC SENSORS



Standard Sensor



High Speed Sensor

## VIDEO OPTIONS




Atomic Force Microscope




Zoom Microscope



# LABORATORY QUALITY RESULTS IN ANY LOCATION



**X - Y SCAN AREA**  
25 x 25mm Motorized



**HEIGHT RANGE**  
2.5mm to 25mm




**DESKTOP DIMENSIONS**  
20 x 30 x 17cm



**SCAN SPEED**  
20mm/s

# PS50 COMPACT OPTICAL PROFILER

- Most advanced compact profilometer
- Small and simple footprint
- Measurement capabilities up to 50mm x 50mm
- All testing capabilities in compact version



## 3D CHROMATIC SENSORS



Standard Sensor



High Speed Sensor

## VIDEO OPTIONS





Atomic Force Microscope




Zoom Microscope

# MOST ADVANCED **COMPACT BENCHTOP**



**X - Y SCAN AREA**  
50 x 50mm Motorized



**HEIGHT RANGE**  
2.5nm to 25mm



**DESKTOP DIMENSIONS**  
38 x 33 x 43cm



**SCAN SPEED**  
20mm/s

# JR100 PORTABLE & HIGH SPEED OPTICAL PROFILER

- Fast measurement (without stitching) using a 100 mm XY axis travel
- Z stage allows setup of measurements at various starting heights
- A high speed sensor gives ultra fast measurements at 382,000 points per second.
- Powerful for quality control



## 3D CHROMATIC SENSORS



Standard Sensor



High Speed Sensor

## VIDEO OPTIONS





Atomic Force Microscope




Zoom Microscope

# PORTABILITY AND HIGH SPEED



**X - Y SCAN AREA**  
100 x 100mm Motorized



**HEIGHT RANGE**  
2.5mm to 25mm



**DESKTOP DIMENSIONS**  
44 x 49 x 32cm



**SCAN SPEED**  
20mm/s

# AFMPRO OPTICAL PROFILER

- 150 x 200mm XY stages and an adjustable height clearance of up to 140mm
- High magnification microscopy
- AFM expands the 3D capabilities into the sub nanometer range
- AFM gives the best lateral accuracy compared to optical techniques
- Easy to select zones on the video to be scanned



## 3D CHROMATIC SENSORS



Standard Sensor



High Speed Sensor

## VIDEO OPTIONS





Atomic Force Microscope




Zoom Microscope

# OPTICAL PROFILER WITH AFM MODULE



**X - Y SCAN AREA**  
200 x 150mm Motorized



**HEIGHT RANGE**  
2.5nm to 25mm




**DESKTOP DIMENSIONS**  
64 x 64 x 82cm



**SCAN SPEED**  
40mm/s

# HS2000 BEST STABILITY OPTICAL PROFILER

- Granite base provides superior stability
- Automated inspection for quality control
- Workstation included to create fully contained stand alone instrument
- Excellent for roughness measurements, combined with advanced automation features



## 3D CHROMATIC SENSORS



Standard Sensor



High Speed Sensor

## VIDEO OPTIONS




Atomic Force Microscope




Zoom Microscope



# HIGH SPEED AND PRECISION FLATNESS TOOL



**X - Y SCAN AREA**  
400 x 500mm Motorized



**HEIGHT RANGE**  
2.5nm to 25mm



**DESKTOP DIMENSIONS**  
101 x 106 x 195cm



**SCAN SPEED**  
500mm/s


***TECHNIQUE***



# HOW IT WORKS

Chromatic Light Technology operates via a process that utilizes white light and a series of spherochromatic lenses. The spherochromatic lenses split the white light into individual wavelengths with unique vertical focal points (vertical distance from surface or height). All wavelengths and their corresponding heights make up the height range measurement scale of a sensor.

The wavelength with the highest intensity will be detected by the spectrometer which processes the wavelength's associated height. During a full raster scan, this process takes a fraction of a second and produces an accurate height map of the surface of interest.



**NO COMPLEX ALGORITHMS** ♦ **NO LEVELING REQUIRED** ♦ **NO X-Y DATA STITCHING**

# THE PROBLEM WITH OTHER TECHNIQUES

## LATERAL RESOLUTION vs LATERAL ACCURACY



**THEM**

**Camera Pixel Size** or **Display Resolution** is often defined as **lateral resolution** to impress clients. Instruments that use camera pixel-based technology require complex algorithms to determine the focal point of the instrument which is problematic for complex surfaces.

**US**

**Chromatic Light** provides **lateral accuracy** which is determined by physics and is directly related to the spot size of the chromatic light source of the optical sensor.

# LASER SCANNING CONFOCAL MICROSCOPE



LASER RADIATION

## **HEALTH HAZARD**

Exposure to laser light reflectivity

## **INCONSISTENT LASER LIGHT WAVELENGTH**

Inconsistencies in wavelength during scanning  
affect accuracy of results

## **DECEPTIVE 'DISPLAY RESOLUTION'**

Lateral & height accuracy are fixed by the objective lens  
making 'Display Resolution' insignificant

## **COMPLEX ALGORITHMS**

Alpha blending algorithms stitch collected data  
layer by layer, grounding accuracy on complex calculations

## **STITCHING REQUIRED**

Objective lenses have limited fixed fields of view  
Stitching of larger areas compromises accuracy of the scan

## **50x SLOWER**

Data acquisition speed up to 7.9 KHz

VS

# CHROMATIC LIGHT OPTICAL SENSOR

## **SAFE WHITE LIGHT**

No need for protective wear

## **UNIFORM & BROAD WHITE LIGHT SPECTRUM**

Changes in wavelength are the data being collected

## **INDEPENDENT LATERAL & HEIGHT ACCURACY**

Lateral & height accuracy can be mixed and matched  
to meet a broad range of scanning requirements

## **NO ALGORITHMS**

Physical wavelength reflected from the surface  
is measured directly for an accurate representative height map

## **NO STITCHING**

Data points are collected continuously providing  
the same level of accuracy for both small and large areas

## **50x FASTER**

Data acquisition speed up to 384 KHz

# LASER MICROSCOPE

## LATERAL ACCURACY

# OPTICAL SENSOR

For 50x objective (370 x 277  $\mu\text{m}$ )

$\pm 2\%$  of measuring value

$\pm 2\% \times 370 \mu\text{m}$

$\approx 15 \mu\text{m}$

w/ stitching algorithms  $\gg 15 \mu\text{m}$



Step size:

$= 5 \mu\text{m}$

**3x BETTER LATERAL ACCURACY**

## HEIGHT ACCURACY

$\approx 0.2 + L/100 \mu\text{m}$

$\approx 0.2 + 950/100 \mu\text{m}$

$\approx 9.7 \mu\text{m}$



950  $\mu\text{m}$  range

$\approx 0.6 \mu\text{m}$

**16x BETTER HEIGHT ACCURACY**

## AREA TESTED

### STITCHING REQUIRED

# scans (25 x 25 mm)

25 000  $\mu\text{m}$  / 370  $\mu\text{m}$  x 25 000  $\mu\text{m}$  / 277  $\mu\text{m}$

68 x 91

$= 6188$  scans



### NO STITCHING

Consistent accuracy across any measurement size

**1 SCAN**


## TEST TIME

6 sec per scan

+ 4 sec displacement & stitching

$= 10 \text{ sec/scan} \times 6188 \text{ scans}$

$= 61860 \text{ seconds}$  ( $\approx 17$  hours)




Scan time (25 x 25 mm)

$= 29.6 \text{ seconds}$



**2090x FASTER**

***TESTING  
SOLUTIONS***



# ROUGHNESS | FINISH

- One second Ra measurement
- Any materials or surface complexity (3D or 2D)
- Automotive roughness finish standards





## ROUGHNESS | FINISH ANALYSIS

- Ra | Sa profile & surface average roughness
- Rq | Sq profile & surface rms roughness
- Rz | Sz maximum height
- Sp | Sv maximum peak & pit height
- SKu | Ssk kurtosis & skewness of height distribution
- Bearing ratio and index
- Sk kernel roughness depth
- Spk | Svk reduced peak height & valley depth
- Sr1 | Sr2 upper & lower material ratio
- Sci & Svi core & valley fluid retention index


AND MORE

# TEXTURE

- Isotropic & anisotropic surfaces
- Hills and valleys analysis



## Texture Direction



## TEXTURE ANALYSIS


- % of isotropy
- 1st, 2nd and 3rd direction
- % of periodicity
- Period
- Density of peaks
- Peak curvature (pointed or rounded)
- Average area of valleys & hills
- Average volume of valleys & hills

AND MORE



# FLATNESS | WARPAGE

- Flatness <1µm over 500mm with no correction




## FLATNESS | WARPAGE ANALYSIS

- 3D & 2D surface waviness & flatness
- Best polynomial match
- Material & bearing ratios
- Distance measurement
- FLTt peak to valley flatness deviation of the surface
- FLTp peak to reference flatness deviation
- FLTv reference to valley flatness deviation
- FLTq rms flatness deviation

AND MORE

# VOLUME | AREA

- Surface subtraction & volume lost
- Corrosion analysis
- Motif and grain analysis





## VOLUME | AREA ANALYSIS

- Volume of void, hills or valleys
- Sdar | Spar developed surface area & projected area
- Volume of void & material from given height
- Map area above or below given heights (%µm²)
- Mean thickness of void & material from given height
- # of grains & average size
- Area & perimeter of grains
- Height, area, volume of motifs
- Max and min pitch of motifs

AND MORE

# GEOMETRY AND SHAPE

- Direct comparison to CAD geometry
- Curvature, radius, angles
- Lateral dimension
- Drill bit studies
- Cutting tools studies





## GEOMETRY AND SHAPE ANALYSIS

- Radius of curvature
- Relative angle measurement
- Distance measurement
- Mean diameter
- Contour analysis
- Rake and wedge angle of drill bit
- K symmetry of cutting edge
- S alpha and gamma dist apex to end of clearance & rake roundness

AND MORE

# STEP HEIGHT | THICKNESS

- Measure through transparent materials
- Transparent film and coating thickness down to 20nm
- Steps from 20nm to 25mm



## STEP HEIGHT | THICKNESS ANALYSIS


- Point to point
- Point to plane
- Maximum, minimum and mean heights
- 3D or 2D map of thickness
- Thickness distribution curve

AND MORE


# MICROSCOPE VIDEO IMAGING

Available on : **ST400, ST500, AFMPRO, & HS2000**

- Ultra zoom lens with coax lighting & detent
- Large area stitching capability
- Color video camera (1200x1600)
- Maximum magnification of 8000X
- Three positions turret (optional)



Broadview map selection tool



PRVision for machine vision capability







Image area selection measurement and image overlay



# ATOMIC FORCE MICROSCOPE

Available on : **AFMPRO**

- Scan of XY 110 $\mu\text{m}$  | high resolution XY 25 $\mu\text{m}$
- Lateral resolution 1.7nm
- Static, dynamic and extended modes
- Max Z range 22 $\mu\text{m}$  | 5 $\mu\text{m}$
- Integrated video camera
- AFM to/from indenter position or video imaging with accuracy of < 0.2 $\mu\text{m}$



# ADVANCED AUTOMATION



- Automatic focus (optical and microscope)
- Automatic analysis template
- Multi sample handling macros
- Easy selection of area under microscope
- Automatic dual frequency for surfaces with varying reflections
- Custom mounting setup of sensors for inline roughness QC
- Pattern recognition, database communications, pass/fail limits



# VISIT OUR APPLICATION NOTES LIBRARY

[nanovea.com/app-notes](https://nanovea.com/app-notes)

Nanovea Optical Profilers **measure any material with a wider range of measurement** than any other Profilometer.



BASE	Jr25	Jr100	PS50	ST400	ST500	HS2000
Type	Portable	Portable & Fast	Compact	Standard	Large Area	Zero Noise / Flatness
X-Y Stage Travel	25 x 25mm	100 x 100mm	50 x 50mm	200 x 150mm	400 x 400mm	400 x 500mm
Z Axis	30mm Manual	25mm Manual	30mm Manual	50mm Motorized	50mm Motorized	100mm Motorized
Maximum X-Y Speed	20 mm/s	20 mm/s	20mm/s	40mm/s	200mm/s	500mm/s
System Dimensions	20 x 30 x 17cm	44 x 49 x 32cm	38 x 33 x 43cm	62 x 62 x 82cm	97 x 72 x 92cm	101 x 106 x 195cm
Rotational Options	N/A	N/A	N/A	Stage or Cylinder	Stage or Cylinder	Software
Video Microscope	N/A	N/A	N/A	Available	Available	Available
Max Sample Weight	No Limit	No Limit	8Kg	8Kg	4-8Kg	7Kg
High Speed Line Sensor	N/A	Included	N/A	Available	Available	Available
Customizable	50mm Stage Travel	N/A	N/A	4 axis & AFM	4 axis	5 axis

### MEASUREMENT TECHNOLOGY

Technique	Non Contact • Chromatic Light
Data Stitching	Not Required within X-Y Stage Travel
Materials Types	ALL - Including Dark, Transparent, & Reflective
Max Surface Angle	Up To 87°

More Information at  
[nanovea.com/profilometers](http://nanovea.com/profilometers)

STANDARD SENSOR (Single Point)	PS1	PS2	PS3	PS4	PS5	PS6
Max Height Range	110µm	300µm	1.1mm	3.5mm	10mm	24mm
Working Distance	3.3mm	10.8mm	12.2mm	16.5mm	26.6mm	20mm
Lateral Accuracy (X-Y)	0.8µm	1.7µm	2.6µm	4.6µm	11.0µm	11.0µm
Height Repeatability (Ra) *	1.9nm	5.4nm	15.8nm	31.6nm	117.0 nm	237.2 nm

HIGH SPEED SENSOR (192 Points)	LS1	LS2	LS3
Max Height Range	200µm	0.95mm	3.9mm
Working Distance	5.3mm	18.5mm	41mm
Height Repeatability (Ra) *	14nm	21nm	70nm
Line Width	0.96mm	1.91mm	4.78mm
Pitch	5µm	10µm	25µm
Lateral Accuracy of each point	1µm	2µm	5µm
Acquisition Rate (points per second)	384KHz	384KHz	384KHz

\* Fixed point on glass, average height variation for 1200points (100 sampling)

Specifications are continuously improving, please contact **NANOVEA** for the latest.

# **N** Today's Standard For Tomorrow's Materials.

**NANOVEA** instruments can be found in renowned education and industrial [organizations around the world](#).

From [aerospace applications](#) to [medical devices](#), thousands of clients at the frontiers of the most demanding industries, with no room for error, rely on our instruments' [unmatched accuracy](#) and [technical superiority](#).



[Nanovea.com](#)

[Nanovea linkedin](#)

[Nanovea Films](#)

[Nanovea FB](#)

[@nanoveainc](#)

